

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Dimethyl sulfoxide - for molecular biology

Product Number: D8418 CAS Number: 67-68-5

MDL: MFCD00002089

Formula: C2H6OS Formula Weight: 78.13 g/mol O H₃C^{ŹŠ}∖CH₃

TEST Specification

Appearance (Color)

Appearance (Form)
Appearance (Clarity)
Proton NMR spectrum
Water (by Karl Fischer)

UV Absorbance Purity (GC)

DNAse, Exonuclease Detection

RNAse Detection Titratable acid

Recommended Retest Period

2 Years

Specification: PRD.2.ZQ5.10000016713

Colorless

Liquid Clear

Conforms to Structure

≤ 0.1 %Pass≥ 99.9 %None DetectedNone Detected

< 0.001 meq

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.